

Product datasheet

Recombinant Human SNRPB2 protein ab107134

画像数 1

製品の詳細

製品名	Recombinant Human SNRPB2 protein
精製度	> 90 % SDS-PAGE. ab107134 is purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	P08579
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHSSGLVPRGSHMDIRPNHTIYINNMDKIKKEE LKRSLYAL FSQFGHVVDIVALKTMKMRGQAFVIFKELGSSTNALRQLQGF PFYGKPMR IQYAKTDSDIISKMRGTFADKEKKKEKKAKTVEQTATTTNK KPGQGTPN SANTQGNSTPNPQVPDYPPNYILFLNNLPEETNEMMLSMLFN QFPGFKEV RLVPRGHDIAFVEFENDGQAGAARDALQGFKITPSHAMKITY AKK
予測される分子量	28 kDa including tags
領域	1 to 225
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab107134** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.316% Tris HCl, 40% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報

機能

Involved in pre-mRNA splicing. This protein is associated with snRNP U2. It binds stem loop IV of U2 snRNA only in presence of the U2A' protein.

配列類似性


Belongs to the RRM U1 A/B" family.

Contains 2 RRM (RNA recognition motif) domains.

細胞内局在

Nucleus.

画像



15% SDS-PAGE showing ab107134 (3µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors