

### Recombinant human SRMS protein ab107697

画像数 2

#### 製品の詳細

製品名	Recombinant human SRMS protein
生理活性	Specific activity: 24 nmol/min/mg
精製度	> 95 % SDS-PAGE. Determined by densitometry. Affinity purified.
発現系	Baculovirus infected Sf9 cells
アクセッション番号	<b><u>Q9H3Y6</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	54 kDa including tags
領域	215 to 488

#### 特性

Our **Abpromise guarantee** covers the use of **ab107697** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Functional Studies
製品の状態	Liquid
備考	<b><u>ab204877</u></b> (Poly (4:1 Glu, Tyr) peptide) can be utilized as a substrate for assessing kinase activity


#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.307% Glutathione, 0.00174% PMSF, 0.00385% DTT, 0.79% Tris HCl, 0.00292% EDTA, 25% Glycerol (glycerin, glycerine), 0.87% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--

## 関連情報


機能	May be involved in proliferation or differentiation of keratinocytes in the skin.
配列類似性	Belongs to the protein kinase superfamily. Tyr protein kinase family. SRC subfamily. Contains 1 protein kinase domain. Contains 1 SH2 domain. Contains 1 SH3 domain.

## 画像



SDS Page Analysis of ab107697

SDS-PAGE - Recombinant human SRMS protein  
(ab107697)



ab107697 specific activity was determined to be 24 nmol/min/mg

Functional Studies - Recombinant human SRMS  
protein (ab107697)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors