

### Recombinant Human Stanniocalcin 2/STC-2 protein ab63281

画像数 1

#### 製品の詳細

製品名	Recombinant Human Stanniocalcin 2/STC-2 protein
精製度	> 90 % SDS-PAGE. Please filter the product by an appropriate sterile filter before using it in cell culture.
エンドキシン・レベル	< 0.100 Eu/μg
発現系	HEK 293 cells
アクセッション番号	<b><u>Q76061</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	TDATNPPEGPQDRSSQQKGRSLSLQNTAEIQHCLVNAGDVGCG VFECFENN SCEIRGLHGICMTFLHNAGKFDAQGKSFIKDALKCKAHALRH RFGCISRK CPAIREMVSQLQRECYLKHDLCAAAQENTRVIVEMIHFKDLL LHEPYVDL VNLLLTGEEVKEAITHSVQVQCEQNWGSLCSILSFCTSAIQ KPPTAPPE RQPQVDRTKLSRAHHGEAGHHLPEPSSRETGRGAKGERGSKS HPNAHARG RVGGLGAQGPSGSSEWEDEQSEYSDIRRAAADYKDDDDK
予測される分子量	33 kDa
領域	25 to 302
配列の追加情報	C-Terminal Flag-tag

#### 特性


Our **Abpromise guarantee** covers the use of **ab63281** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	ELISA

製品の状態	Lyophilized
備考	Previously labelled as Stanniocalcin 2.
<b>前処理および保存</b>	
保存方法および安定性	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituents: 0.242% Tris, 0.116% Sodium chloride
再構成	Add deionized water to prepare a working stock solution of approximately 0.5 mg/mL and let the lyophilized pellet dissolve completely.
<b>関連情報</b>	
機能	Has an anti-hypocalcemic action on calcium and phosphate homeostasis.
組織特異性	Expressed in a variety of tissues including muscle, heart, pancreas, kidney, spleen, prostate, small intestine, colon and peripheral blood leukocytes.
配列類似性	Belongs to the stanniocalcin family.
細胞内局在	Secreted.

## 画像



LANE 1. M.W. marker - 14, 21, 31, 45, 66, 97 kDa. LANE 2. reduced and heated sample, 5µg/lane. LANE 3. non-reduced and non-heated sample, 5µg/lane.

SDS-PAGE - Recombinant Stanniocalcin 2/STC-2 protein (ab63281)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors