

Recombinant Human TIA1 protein ab131879

画像数 1

製品の詳細

製品名	Recombinant Human TIA1 protein
発現系	Wheat germ
アクセッション番号	<u>Q96B58</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MEDEMPKTL YVGNLSRDVTEALILQLFSQIGPCKNCKMIMDT AGNDPYCF VEFHEHRHAAAALAAMNGRKIMGKEVKVNWATTPSSQKKDTS SSTVVSTQ RSQDHFHVFGDLSPEITTEDIKAAFAPFGRISDARVVKDMA TGKSKGYG FVSFFNKWDAENAIQQMGQWLGGQRQIRTNWATRKPPAPKST YECRCIGE EKEMWNFGKEYARF
予測される分子量	49 kDa including tags
領域	1 to 214

特性

Our **Abpromise guarantee** covers the use of **ab131879** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE
	ELISA
製品の状態	Liquid

前処理および保存


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
------------	---

Constituents: 0.31% Glutathione, 0.79% Tris HCl

関連情報

機能	Involved in alternative pre-RNA splicing and regulation of mRNA translation by binding to AU-rich elements (AREs) located in mRNA 3' untranslated regions (3' UTRs). Possesses nucleolytic activity against cytotoxic lymphocyte target cells. May be involved in apoptosis.
配列類似性	Contains 3 RRM (RNA recognition motif) domains.
細胞内局在	Cytoplasmic granule. Nucleus. Accumulates in cytoplasmic stress granules (SG) following cellular damage.

画像



12.5% SDS-PAGE stained with Coomassie Blue showing ab131879 at approximately 49.28 kDa.

SDS-PAGE - Recombinant Human TIA1 protein
(ab131879)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors