

## Product datasheet

# Recombinant Mouse 4-1BBL protein (His tag) ab271407

画像数 2

### 製品の詳細

---

製品名	Recombinant Mouse 4-1BBL protein (His tag)
精製度	>= 90 % SDS-PAGE.
発現系	HEK 293 cells
アクセッション番号	<b>P41274</b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Mouse
配列	RTEPRPALT I TTSPNLGTRE NNADQVTPVS HIGCPNTTQQ GSPVFAKLLA KNQASLCNTT LNWHSQDGAG SSYLSQGLRY EEDKKELVVD SPGLYYVFLE LKLSPTFTNT GHKVQGWVSL VLQAKPQVDD FDNLALTVEL FPCSMENKLV DRSWSQLLLL KAGHRLSVGL RAYLHGAQDA YRDWELSYPN TTSFGLFLVK PDNPWE
予測される分子量	26 kDa including tags
分子量情報	Protein MW is higher than predicted on SDS-PAGE due to glycosylation.
領域	104 to 309
タグ	His tag N-Terminus
配列の追加情報	Genbank: NM_009404

### 特性

---

Our **Abpromise guarantee** covers the use of **ab271407** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE

製品の状態	Liquid
-------	--------

前処理および保存

**保存方法および安定性**

Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituents: 0.13% Sodium phosphate, 0.64% Sodium chloride, 0.02% Potassium chloride, 20% Glycerol (glycerin, glycerine)

**関連情報**

**機能**

Cytokine that binds to TNFRSF9. Induces the proliferation of activated peripheral blood T-cells. May have a role in activation-induced cell death (AICD). May play a role in cognate interactions between T-cells and B-cells/macrophages.

**組織特異性**

Expressed in brain, placenta, lung, skeletal muscle and kidney.


**配列類似性**

Belongs to the tumor necrosis factor family.

**細胞内局在**


Membrane.

**画像**



SDS-PAGE analysis of ab271407.

SDS-PAGE - Recombinant Mouse 4-1BBL protein (Tagged) (ab271407)



All lanes :

Lane 1 : ab271407

Lane 2 : MW markers

Predicted band size: 27 kDa

Western blot analysis of ab271407.

Western blot - Recombinant Mouse 4-1BBL protein (His tag) (ab271407)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors