

# Recombinant mouse Alpha-synuclein protein aggregate Type 1 (Active) ab246002

4 References [画像数 6](#)

### 製品の詳細

製品名	Recombinant mouse Alpha-synuclein protein aggregate Type 1 (Active)
生理活性	Endogenous alpha-synuclein phosphorylation. 100 µM alpha synuclein protein monomer seeded with 10 nM alpha synuclein protein PFF (ab246002) in 25 µM Thioflavin T (PBS pH 7.4, 100 µl reaction volume) generated an increased fluorescence intensity after incubation at 37°C with shaking at 600 rpm for 24 hours. Fluorescence was measured by excitation at 450 nm and emission at 485 nm on a microplate reader.
精製度	> 95 % SDS-PAGE. Ion-exchange purified.
エンドトキシン・レベル	< 5.000 Eu/ml
発現系	Escherichia coli
アクセッション番号	<b><u>O55042</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Mouse
配列	MDVFMKGLSKAKEGVVAAAEKTKQGVAAEAGKTKEGVLYVGS KTKEGVVH GVTTVAEKTKEQVTNVGGAVVTGVTAVAQKTVEGAGNIAAAT GFVKKDQM GKGEEGYPQEGILEDMPVDPGSEAYEMPSEEGYQDYEPEA
予測される分子量	15 kDa
領域	1 to 140
配列の追加情報	NP_001035916.1


### 特性

Our **Abpromise guarantee** covers the use of **ab246002** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Immunocytochemistry
	Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections)

	Functional Studies
	Electron Microscopy
	SDS-PAGE
製品の状態	Liquid
備考	Active Mouse recombinant Alpha Synuclein Pre-Formed Fibrils (Type 1). For best results, sonicate immediately prior to use.
<hr/>	
前処理および保存	
<hr/>	
保存方法および安定性	Shipped on Dry Ice. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituent: 95% PBS This product is an active protein and may elicit a biological response in vivo, handle with caution.
<hr/>	
関連情報	
<hr/>	
機能	May be involved in the regulation of dopamine release and transport. Induces fibrillization of microtubule-associated protein tau. Reduces neuronal responsiveness to various apoptotic stimuli, leading to a decreased caspase-3 activation.
組織特異性	Expressed principally in brain but is also expressed in low concentrations in all tissues examined except in liver. Concentrated in presynaptic nerve terminals.
関連疾患	Genetic alterations of SNCA resulting in aberrant polymerization into fibrils, are associated with several neurodegenerative diseases (synucleinopathies). SNCA fibrillar aggregates represent the major non A-beta component of Alzheimer disease amyloid plaque, and a major component of Lewy body inclusions. They are also found within Lewy body (LB)-like intraneuronal inclusions, glial inclusions and axonal spheroids in neurodegeneration with brain iron accumulation type 1. Parkinson disease 1 Parkinson disease 4 Dementia Lewy body
配列類似性	Belongs to the synuclein family.
ドメイン	The 'non A-beta component of Alzheimer disease amyloid plaque' domain (NAC domain) is involved in fibrils formation. The middle hydrophobic region forms the core of the filaments. The C-terminus may regulate aggregation and determine the diameter of the filaments.
翻訳後修飾	Phosphorylated, predominantly on serine residues. Phosphorylation by CK1 appears to occur on residues distinct from the residue phosphorylated by other kinases. Phosphorylation of Ser-129 is selective and extensive in synucleinopathy lesions. In vitro, phosphorylation at Ser-129 promoted insoluble fibril formation. Phosphorylated on Tyr-125 by a PTK2B-dependent pathway upon osmotic stress. Hallmark lesions of neurodegenerative synucleinopathies contain alpha-synuclein that is modified by nitration of tyrosine residues and possibly by dityrosine cross-linking to generated stable oligomers. Ubiquitinated. The predominant conjugate is the diubiquitinated form. Acetylation at Met-1 seems to be important for proper folding and native oligomeric structure.
細胞内局在	Cytoplasm, cytosol. Membrane. Nucleus. Cell junction, synapse. Secreted. Membrane-bound in dopaminergic neurons.




Immunocytochemistry - Recombinant mouse Alpha-synuclein protein aggregate (Active) (ab246002)

Primary rat hippocampal neurons (DV16) show lewy body inclusion formation and loss of cells when treated with ab246002 at 4  $\mu$ g/ml (D-F) on DVI2, but not when treated with a control (A-C). Tissue: Primary hippocampal neurons. Species: Sprague-Dawley rat. Fixation: 3% formaldehyde from PFA for 20 min. Blocker: 1:1 PBS:proprietary block and 30 mL/mL of 0.1% triton-X 100 for 30 min. Primary Antibody: Mouse anti-pSer129 Antibody (1/1000) and Rabbit anti-pSer129 (1/800) for 24 hours at 4°C. Secondary Antibody: ATTO 546 Donkey Anti-Mouse (1/700) and ATTO 488 Donkey Anti-Rabbit (1/700) for 1 hour at room temperature (composite green). Counterstain: Hoechst (blue) nuclear stain at 1/3000 for 1 hour at room temperature. Localization: Lewy body inclusions. Magnification: 20x.




Immunohistochemistry (PFA fixed) - Recombinant mouse Alpha-synuclein protein (Active) (ab246002)

Immunohistochemistry analysis of rat brain injected with ab246002. Species: Female Sprague-Dawley Rat. Rat was injected with 2 $\mu$ L ab246002 in each of 2 injection sites: AP+1.6, ML+2.4, DV-4.2 from skull; and AP-1.4, ML+0.2, DV-2.8 from skull. 30 days post-injection. Fixation: Saline perfusion followed by 4% PFA fixation for 48 hours. Primary antibody: rabbit monoclonal anti-pSer129 alpha synuclein. Secondary Antibody: Biotin-SP Donkey Anti-Rabbit IgG (H+L) at 1/500 for 2 hours in cold room with shaking. ABC signal amplification, DAB staining. Magnification: 20x. Alpha synuclein pathology is seen in the periform/insular cortex and the cingulate cortex on both the same (ipsi) and opposite (contra) sides as the injection sites.




Electron Microscopy - Recombinant mouse Alpha-synuclein protein (Active) (ab246002)

TEM of ab246002. Fibrils were sonicated and image was taken at 100kx magnification.



Electron Microscopy - Recombinant mouse Alpha-synuclein protein (Active) (ab246002)


TEM of ab246002. Image was taken at 100kx magnification.



Functional Studies - Recombinant mouse Alpha-synuclein protein (Active) (ab246002)

ab246002 seeds the formation of new Alpha synuclein fibrils from the pool of active alpha synuclein monomers. Thioflavin T is a fluorescent dye that binds to beta sheet-rich structures, such as those in alpha synuclein fibrils. Upon binding, the emission spectrum of the dye experiences a red-shift, and increased fluorescence intensity. Thioflavin T emission curves show increased fluorescence (correlated to alpha synuclein protein aggregation) over time when 10 nM of ab246002 is combined with 100 µM of active Alpha synuclein monomer, as compared to ab246002 and active alpha Synuclein monomer alone. Thioflavin T ex = 450 nm, em = 485 nm.

SDS-PAGE analysis of ab246002 (2 µg).



SDS-PAGE - Recombinant mouse Alpha-synuclein protein (Active) (ab246002)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors