

### Recombinant Mouse DR6 protein ab109063

#### 製品の詳細

製品名	Recombinant Mouse DR6 protein
精製度	> 95 % SDS-PAGE.
エンドトキシン・レベル	< 0.100 Eu/μg
発現系	HEK 293 cells
アクセッション番号	<b><u>Q9EPU5</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Mouse
予測される分子量	80 kDa
領域	42 to 348

#### 特性

Our **Abpromise guarantee** covers the use of **ab109063** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized
備考	PBS containing at least 0.1% BSA should be used for further dilutions.

#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C. Constituent: PBS
再構成	Reconstitute with 50 μl sterile water.

#### 関連情報

機能	May activate NF-kappa-B and promote apoptosis. May activate JNK and be involved in T-cell
----	---

differentiation. Required for both normal cell body death and axonal pruning. Trophic-factor deprivation triggers the cleavage of surface APP by beta-secretase to release sAPP-beta which is further cleaved to release an N-terminal fragment of APP (N-APP). N-APP binds TNFRSF21 triggering caspase activation and degeneration of both neuronal cell bodies (via caspase-3) and axons (via caspase-6).

**組織特異性**

Highly expressed in heart, brain, placenta, pancreas, lymph node, thymus and prostate. Detected at lower levels in lung, skeletal muscle, kidney, testis, uterus, small intestine, colon, spleen, bone marrow and fetal liver. Very low levels were found in adult liver and peripheral blood leukocytes.

**配列類似性**

Contains 1 death domain.  
Contains 4 TNFR-Cys repeats.

**細胞内局在**

Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors