

Product datasheet

Recombinant Mouse Eotaxin 2 protein ab50164

製品の詳細

製品名	Recombinant Mouse Eotaxin 2 protein
精製度	> 95 % SDS-PAGE. ab50164 purity was assessed also by HPLC. Endotoxin level is less than 0.1 ng per µg(1EU/µg).
エンドトキシン・レベル	< 0.100 EU/µg
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Mouse
配列	VTIPSSCCTS FISKKIPENR VVSYQLANGS ICPKAGVIFI TKKGHKICTD PKLLWVQRHI QKLDKKNQP SKGAKAVRTK FAVQRRRGNS TEV
領域	27 to 119

特性

Our **Abpromise guarantee** covers the use of **ab50164** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Functional Studies
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term.
再構成	Reconstitute to 1mg/ml in distilled water.

関連情報

機能	Chemotactic for resting T-lymphocytes, and eosinophils. Has lower chemotactic activity for
----	--

neutrophils but none for monocytes and activated lymphocytes. Is a strong suppressor of colony formation by a multipotential hematopoietic progenitor cell line. Binds to CCR3.

組織特異性

Activated monocytes and activated T lymphocytes.

配列類似性

Belongs to the intercrine beta (chemokine CC) family.

翻訳後修飾

N-glycosylated.

細胞内局在

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors