

# Recombinant mouse Eph receptor A7/EPHA7 protein (Active) ab268530

画像数 2

### 製品の詳細

製品名	Recombinant mouse Eph receptor A7/EPHA7 protein (Active)
生理活性	The specific activity of ab268530 was 35 nmol/min/mg in a kinase assay using Poly E4Y peptide as substrate.
精製度	> 95 % SDS-PAGE. Affinity purified.
発現系	Baculovirus infected Sf9 cells
アクセッション番号	<b><u>Q8R381</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Mouse
領域	580 to 994
タグ	GST tag N-Terminus
配列の追加情報	BC026153

### 特性

Our **Abpromise guarantee** covers the use of **ab268530** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Functional Studies SDS-PAGE
製品の状態	Liquid

### 前処理および保存


保存方法および安定性	Shipped on Dry Ice. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. pH: 7.50 Constituents: 0.79% Tris HCl, 0.87% Sodium chloride, 0.31% Glutathione, 0.003% EDTA, 0.004% DTT, 0.002% PMSF, 25% Glycerol (glycerin, glycerine)
------------	---

This product is an active protein and may elicit a biological response in vivo, handle with caution.

## 関連情報


機能	Receptor for members of the ephrin-A family. Binds to ephrin-A1, -A2, -A3, -A4 and -A5.
組織特異性	Widely expressed.
配列類似性	Belongs to the protein kinase superfamily. Tyr protein kinase family. Ephrin receptor subfamily. Contains 2 fibronectin type-III domains. Contains 1 protein kinase domain. Contains 1 SAM (sterile alpha motif) domain.
細胞内局在	Membrane.

## 画像



The specific activity of ab268530 was 35 nmol/min/mg in a kinase assay using Poly E4Y peptide as substrate.

Functional Studies - Recombinant mouse Eph receptor A7/EPHA7 protein (Active) (ab268530)



SDS-PAGE analysis of ab268530.

MW ~73 kDa.

SDS-PAGE - Recombinant mouse Eph receptor A7/EPHA7 protein (Active) (ab268530)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors