

Product datasheet

Recombinant mouse Rab3A protein ab90790

製品の詳細

製品名	Recombinant mouse Rab3A protein
生理活性	Activity: 100 pmol of protein can bind > 80 pmol of GDP.
精製度	> 90 % SDS-PAGE.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Mouse

特性

Our **Abpromise guarantee** covers the use of **ab90790** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

製品の状態 Liquid

備考

The proprietary tag facilitates the protein's application in typical proprietary tag pull down assays.

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.095% Magnesium chloride, 0.0156% Beta mercaptoethanol, 0.79% Tris HCl, 0.0292% EDTA, 0.58% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	---

関連情報

機能	Involved in exocytosis by regulating a late step in synaptic vesicle fusion. Could play a role in
----	---

組織特異性

neurotransmitter release by regulating membrane flow in the nerve terminal.

配列類似性

Specifically expressed in brain.

細胞内局在

Belongs to the small GTPase superfamily. Rab family.

Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors