

Recombinant OV-16 antigen protein (His tag) ab272024

画像数 1

製品の詳細

| | |
|-------------|---|
| 製品名 | Recombinant OV-16 antigen protein (His tag) |
| 精製度 | > 90 % SDS-PAGE. |
| 発現系 | Escherichia coli |
| タンパク質長 | Full length protein |
| Animal free | No |
| 由来 | Recombinant |

特性


| | |
|--|-------------|
| Our Abpromise guarantee covers the use of ab272024 in the following tested applications. | |
| The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user. | |
| アプリケーション | SDS-PAGE |
| 製品の状態 | Lyophilized |

前処理および保存

| | |
|------------|---|
| 保存方法および安定性 | Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: Tris buffer, 6% Trehalose |
|------------|---|

関連情報

画像



SDS-PAGE - Recombinant OV-16 antigen protein
(His tag) (ab272024)

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) analysis
of ab272024 with 5% enrichment gel and 15% separation gel.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors