

Recombinant rat IL-17A protein ab78597

画像数 2

製品の詳細

製品名	Recombinant rat IL-17A protein
生理活性	Biological Activity: ab78597 is fully biologically active when compared to standards via dose-dependent induction of IL6 production in cultured mouse NIH 3T3 fibroblasts. The ED ₅₀ for ab78597 in this assay was 30ng/ml.
精製度	> 95 % SDS-PAGE. Purity is greater than 97% as determined by reducing and non-reducing SDS-PAGE. Endotoxin level is <0.01ng/ug or <0.1EU/ug. This product is filter sterilised.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Rat
配列	MAVLIPQSSV CPNAEANNFL QNVKVNLIKVL NSLSSKASSR RPSDYLNIRST SPWTLSRNED PDRYPSVIWE AQCRHQRCVN AEGKLDHMMN SVLIQQEILV LKREPEKCPF TFRVEKMLVG VGCTCVSSIV RHAS

特性

Our **Abpromise guarantee** covers the use of **ab78597** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized

前処理および保存


保存方法および安定性	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C. Constituent: 0.294% Sodium citrate This product is an active protein and may elicit a biological response in vivo, handle with caution.
再構成	Reconstitute with sterile water at a concentration of 0.1 mg/ml, which can be further diluted into

other aqueous solutions. Reconstituted material should be aliquoted and frozen at -20C. It is recommended to add a carrier protein (0.1% HSA or BSA) for long term storage.


関連情報

機能	Induces stromal cells to produce proinflammatory and hematopoietic cytokines. Enhances the surface expression of the intracellular adhesion molecule-1 (ICAM-1) in fibroblasts.
組織特異性	Restricted to activated memory T-cells.
配列類似性	Belongs to the IL-17 family.
翻訳後修飾	Found both in glycosylated and nonglycosylated forms.
細胞内局在	Secreted.

画像



SDS PAGE analysis of ab78597 under non-reducing (-) and reducing (+) conditions. Stained with Coomassie Blue.



ab78597 used in Functional Studies.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

-
- Replacement or refund for products not performing as stated on the datasheet
 - Valid for 12 months from date of delivery
 - Response to your inquiry within 24 hours

 - We provide support in Chinese, English, French, German, Japanese and Spanish
 - Extensive multi-media technical resources to help you
 - We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors