

### Recombinant rat IP10 protein ab75477

画像数 1

#### 製品の詳細

製品名	Recombinant rat IP10 protein
精製度	> 95 % SDS-PAGE. Purity : Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Rat
配列	IPLARTVRCT CIDFHEQPLR PRAIGKLEII PASLSCPHVE IIATMKKNNE KRCLNPESEA IKSLLKAVSQ RRSKRAP
領域	22 to 98

#### 特性

Our **Abpromise guarantee** covers the use of **ab75477** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE
	Functional Studies
製品の状態	Lyophilized


#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle.  This product is an active protein and may elicit a biological response in vivo, handle with caution.
再構成	Reconstitute in water to a concentration of 0.1-1.0 mg/ml. Do not vortex. This solution can be stored at 2-8oC for up to 1 week. For extended storage, it is recommended to further dilute in a buffer containing a carrier protein (example 0.1% BSA) and store in working aliquots at -20oC to -80oC.

## 関連情報

機能	Chemotactic for monocytes and T-lymphocytes. Binds to CXCR3.
配列類似性	Belongs to the intercrine alpha (chemokine CxC) family.
翻訳後修飾	CXCL10(1-73) is produced by proteolytic cleavage after secretion from keratinocytes.
細胞内局在	Secreted.

## 画像



Anti-IP10 antibody (**ab9807**) at 1 µg/ml + Recombinant rat IP10 protein (ab75477) at 0.1 µg

### Secondary

Goat Anti-Rabbit IgG H&L (HRP) preadsorbed (**ab97080**) at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

**Exposure time:** 30 seconds

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors