

### Density Gradient Medium (Iodixanol) ab286850

画像数 1


#### 製品の概要

製品名	Density Gradient Medium (Iodixanol)
特記事項	This product is manufactured by BioVision, an Abcam company and was previously called M1248 OptiPrep™ Density Gradient Medium. M1248-250 is the same size as the 250 ml size of ab286850.  A ready-made sterile solution of 60% iodixanol in water with a density of 1.32 g/ml.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at Room Temperature. Store at Room Temperature. Do Not Freeze. This product is air and light sensitive and impurities can occur as a result of air oxidation or due to metabolism by microbes.
バッファー	Constituent: 60% Iodixanol

#### 画像



Chemical Structure - Density Gradient Medium  
(Iodixanol) (ab286850)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors