

Product datasheet

Normal Rabbit Serum ab166640

8 References

製品の概要

製品名	Normal Rabbit Serum
由来種	Rabbit
特記事項	Rabbit serum for blocking in antibody-based applications or for use in cell culture media. North America origin. Mixed sex rabbits. New Zealand White strain.

製品の特性

製品の状態	Liquid
保存方法	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
関連性	Serum is the clear, straw-colored, liquid portion of blood plasma that does not contain fibrinogen or blood cells and remains fluid after clotting.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors