

Donkey Anti-Mouse IgG H&L (6nm Gold) preadsorbed ab105276

1 References

製品の概要

製品名	Donkey Anti-Mouse IgG H&L (6nm Gold) preadsorbed
由来種	Donkey
ターゲット生物種	Mouse
特異性	Recognizes the heavy chains on mouse IgG and the light chains common to most murine immunoglobulins. Does not crossreact with non-immunoglobulin serum proteins. Donkey Anti-Mouse IgG H&L (6nm Gold) preadsorbed exhibits minimal cross reactivity to Bovine, chicken, goat, guinea pig, hamster, horse, human, rabbit and sheep serum proteins.
アプリケーション	適用あり: Electron Microscopy, WB, ELISA
吸着処理血清	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Sheep more details
標識	Gold 6nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 8.50 Preservative: 0.04% Sodium azide Constituents: 1.34% BSA, 0.25% Sodium borate, 0.8% Sodium chloride, 0.09% Dibasic monohydrogen sodium phosphate, 89.08% Water
精製度	Immunogen affinity purified
特記事項 (精製)	Purified by affinity chromatography.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab105276の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Electron Microscopy		1/20 - 1/40.
WB		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors