

Donkey Anti-Rabbit IgG H&L (Texas Red ®) ab6800

22 References

製品の概要

製品名	Donkey Anti-Rabbit IgG H&L (Texas Red ®)
由来種	Donkey
ターゲット生物種	Rabbit
アプリケーション	適用あり: Immunomicroscopy, Flow Cyt, IHC-P, IHC-Fr, ICC/IF, ELISA
免疫原	Rabbit IgG, whole molecule
標識	Texas Red ®. Ex: 596nm, Em: 620nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 6.50 Preservative: 0.01% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Tripotassium orthophosphate
精製度	Affinity purified
特記事項 (精製)	Texas Red Sulfonyl Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm Emission Wavelength: 620 nm Fluorochrome/Protein Ratio: 1.7 moles Texas Red per mole of Goat IgG
特記事項 (標識)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose
ポリ/モノ	ポリクローナル

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab6800の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Immunomicroscopy		Use at an assay dependent dilution.
Flow Cyt		Use at an assay dependent dilution.
IHC-P		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.
ICC/IF		1/1000 - 1/5000.
ELISA		1/10000 - 1/50000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors