

Donkey Anti-Rat IgG H&L (HRP) ab102182

7 References

製品の概要

| | |
|----------|---|
| 製品名 | Donkey Anti-Rat IgG H&L (HRP) |
| 由来種 | Donkey |
| ターゲット生物種 | Rat |
| 特異性 | By immunoelectrophoresis and ELISA this antibody reacts specifically with rat IgG and with light chains common to other rat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with IgG from other species. |
| アプリケーション | 適用あり: WB, ELISA, ICC, IHC-P |
| 標識 | HRP |

製品の特性

| | |
|-----------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative |
| 精製度 | Immunogen affinity purified |
| 特記事項 (精製) | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Horse Radish Peroxidase. |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab102182の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| WB | | 1/5000 - 1/50000. (Colorimetric: 1/5000 - 1/30000; Chemiluminescence: 1/10000 - 1/50000). |

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--------------------------------|
| ELISA | | 1/10000 - 1/100000. (Primary). |
| ICC | | 1/200 - 1/5000. |
| IHC-P | | 1/200 - 1/5000. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors