

Goat Anti-Chicken IgY H&L (Alkaline Phosphatase) preadsorbed  
ab7119

製品の概要

製品名	Goat Anti-Chicken IgY H&L (Alkaline Phosphatase) preadsorbed	
由来種	Goat	
ターゲット生物種	Chicken	
アプリケーション	適用あり: ICC/IF, Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, WB	
吸着処理血清	Cow, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Rat, Sheep	<a href="#">more details</a>
免疫原	Full length native chicken IgY (purified).	
標識	Alkaline Phosphatase	

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Constituents: 0.00136% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride
精製度	Affinity purified
特記事項 (精製)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using chicken IgG coupled to agarose beads.
特記事項 (標識)	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
ポリ/モノ	ポリクローナル

アプリケーション

The Abpromise guarantee      Abpromise保証は、 次のテスト済みアプリケーションにおけるab7119の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		

アプリケーション	Abreviews	特記事項
Dot blot		
ELISA		
IHC-P		
IHC-Fr		
Immunomicroscopy		
WB		

#### 追加情報

Dot: Use at an assay dependant dilution.

ELISA: 1/1000 - 1/4000.

IHC-P: Use at an assay dependant dilution.

IHC-Fr: Use at an assay dependant dilution.

IM: Use at an assay dependant dilution.

WB: Use at an assay dependant dilution.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors