

# Goat Anti-Chicken IgY H&L (DyLight® 488) preadsorbed ab96951

**6 References** [画像数 1](#)

### 製品の概要

製品名	Goat Anti-Chicken IgY H&L (DyLight® 488) preadsorbed
由来種	Goat
ターゲット生物種	Chicken
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Chicken IgG and with light chains common to other Chicken immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. Reduced cross-reactivity to bovine, goat, horse, human, mouse, pig, rabbit and rat was detected.
アプリケーション	<b>適用あり:</b> IHC-P, ICC/IF, Flow Cyt
吸着処理血清	Cow, Goat, Horse, Human, Pig, Rabbit, Rat <a href="#">more details</a>
標識	DyLight® 488. Ex: 493nm, Em: 518nm

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS
精製度	Immunogen affinity purified
特記事項(精製)	Antiserum was cross adsorbed using bovine, horse, human, mouse, pig, rabbit and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 488.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション


## The Abpromise guarantee

**Abpromise保証は、次のテスト済みアプリケーションにおけるab96951の使用に適用されます**

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/500.
ICC/IF		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

## 画像



Goat Anti-Chicken IgY H&L (DyLight® 488)  
preadsorbed (ab96951)

Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors