

Goat F(ab) Anti-Mouse IgG H&L ab6668

36 References 画像数 2

製品の概要

製品名	Goat F(ab) Anti-Mouse IgG H&L
由来種	Goat
ターゲット生物種	Mouse
アプリケーション	適用あり: WB, ICC/IF, IHC, ELISA
免疫原	Mouse IgG whole molecule

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride
特記事項 (精製)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, papain digestion and chromatographic separation. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum. No reaction was observed against anti-Papain or anti-Goat IgG F(c).”
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab6668の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/2000 - 1/10000.
ICC/IF		1/1000 - 1/5000.

アプリケーション	Abreviews	特記事項
IHC		1/1000 - 1/5000.
ELISA		1/50000.

画像



SDS-PAGE results of Goat Fab Anti-Mouse IgG (H&L) ab6668.


Lane 1: Non-reduced Goat Fab Anti-Mouse.

Lane 2: Opal Prestained Ladder.

Lane 3: Reduced Goat Fab Anti-Mouse.

Load: 1.0 µg.

Coomassie stained and imaged.



Dot Blot results of Goat Fab Anti-Mouse IgG Antibody ab6668.

Dots are Mouse IgG at (1) 100ng, (2) 33.3ng, (3) 11.1ng, (4) 3.70ng, (5) 1.23ng.

Blocked for 30 min at RT.

Primary Antibody: Goat Fab Anti-Mouse IgG at 1 µg/mL for 1 hr at RT.

Secondary Antibody: Donkey Anti-Goat IgG HRP at 1:40,000 for 30 min at RT.

Imaged with BioRad ChemiDoc, Chemi filter.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors