

## Product datasheet

# Goat F(ab')<sub>2</sub> Anti-Mouse IgG - (Fab')<sub>2</sub> (DyLight® 488), pre-adsorbed ab98757

### 1 References [画像数 1](#)

#### 製品の概要

|          |  |
|----------|--|
| 製品名      | Goat F(ab') <sub>2</sub> Anti-Mouse IgG - (Fab') <sub>2</sub> (DyLight® 488), pre-adsorbed |
| 由来種      | Goat   |
| ターゲット生物種 | Mouse  |
| アプリケーション | <b>適用あり:</b> Flow Cyt, ICC/IF, IHC-P   |
| 吸着処理血清   | Human, Rat <a href="#">more details</a>  |
| 標識       | DyLight® 488. Ex: 493nm, Em: 518nm   |

#### 製品の特性


|          |  |
|----------|--|
| 製品の状態    | Liquid   |
| 保存方法     | Shipped at 4°C. Store at +4°C.   |
| バッファー    | pH: 6.8<br>Preservative: 0.09% Sodium azide<br>Constituents: PBS, 0.2% BSA   |
| 精製度      | Immunogen affinity purified  |
| 特記事項(精製) | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') <sub>2</sub> fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to DyLight® 488. |
| ポリ/モノ    | ポリクローナル  |
| アイソタイプ   | IgG  |

#### アプリケーション

**The Abpromise guarantee** **Abpromise保証は、次のテスト済みアプリケーションにおけるab98757の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項          |
|----------|-----------|---------------|
| Flow Cyt |           | 1/50 - 1/200. |
| ICC/IF   |           | 1/50 - 1/500. |
| IHC-P    |           | 1/50 - 1/500. |

## 画像



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

Goat F(ab')<sub>2</sub> Anti-Mouse IgG - (Fab')<sub>2</sub> (DyLight® 488), pre-adsorbed (ab98757)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors