

Product datasheet

Goat F(ab')₂ Anti-Mouse IgG - (Fab')₂ (DyLight® 550), pre-adsorbed ab98758

[3 References](#) [画像数 1](#)

製品の概要

製品名	Goat F(ab') ₂ Anti-Mouse IgG - (Fab') ₂ (DyLight® 550), pre-adsorbed
由来種	Goat
ターゲット生物種	Mouse
アプリケーション	適用あり: Flow Cyt, ICC/IF, IHC-P
吸着処理血清	Human, Rat more details
標識	DyLight® 550. Ex: 562nm, Em: 576nm

製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') ₂ fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to DyLight® 550.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	DyLight® 550 replaces DyLight® 549.

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab98758の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		1/50 - 1/200.
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.

画像



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors