

Goat F(ab')₂ Anti-Mouse IgG - Fc (AP), pre-adsorbed ab5880

1 References

製品の概要

| | |
|----------|--|
| 製品名 | Goat F(ab') ₂ Anti-Mouse IgG - Fc (AP), pre-adsorbed |
| 由来種 | Goat |
| ターゲット生物種 | Mouse |
| 特異性 | Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Alkaline Phosphatase, anti-Goat Serum, Mouse IgG, Mouse IgG F(c) and Mouse Serum. No reaction was observed against anti-Pepsin, anti-Goat IgG F(c), Mouse IgG F(ab') ₂ or Bovine, Horse and Human Serum Proteins. |
| アプリケーション | 適用あり: ICC/IF, Dot blot, WB, ELISA, IHC |
| 吸着処理血清 | Cow, Horse, Human more details |
| 免疫原 | Mouse IgG F(c) fragment |
| 標識 | Alkaline Phosphatase |

製品の特性

| | |
|-----------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | Constituents: 0.00136% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride BSA is IgG and protease free |
| 特記事項 (精製) | This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation. |
| 特記事項 (標識) | Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons). |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab5880の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|-------------------------------------|
| ICC/IF | | 1/200 - 1/1000. |
| Dot blot | | Use at an assay dependent dilution. |
| WB | | 1/500 - 1/2500. |
| ELISA | | 1/2000 - 1/10000. |
| IHC | | Use at an assay dependent dilution. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors