

Goat F(ab')₂ Anti-Mouse IgG - H&L (AP), pre-adsorbed ab98724

1 References

製品の概要

製品名	Goat F(ab') ₂ Anti-Mouse IgG - H&L (AP), pre-adsorbed
由来種	Goat
ターゲット生物種	Mouse
アプリケーション	適用あり: ICC, IHC-P, ELISA, WB
吸着処理血清	Human, Rat more details
免疫原	Other Immunogen Type corresponding to IgG.
標識	Alkaline Phosphatase

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.58% Sodium chloride, 0.0095% Magnesium chloride, 0.00136% Zinc chloride, 1.19% HEPES, 0.2% BSA
精製度	Immunogen affinity purified
特記事項(精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. F(ab') ₂ fragment were generated using a pepsin digestion. Fc fragments and whole IgG molecules have been removed. Fragments were conjugated to Alkaline Phosphatase.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab98724の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC		Use at an assay dependent dilution.
IHC-P		1/100 - 1/1000.
ELISA		1/1000 - 1/10000. (Primary).
WB		1/1000 - 1/10000. (Chemiluminescent detection).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors