

Goat Anti-Guinea pig IgG H&L (HRP) ab97155

10 References

製品の概要

製品名	Goat Anti-Guinea pig IgG H&L (HRP)
由来種	Goat
ターゲット生物種	Guinea pig
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Guinea Pig IgG and with light chains common to other Guinea Pig immunoglobulins. No antibody was detected against non immunoglobulin serum proteins.
アプリケーション	適用あり: ICC, IHC-P, ELISA, WB
免疫原	Other Immunogen Type corresponding to Guinea pig IgG.
標識	HRP

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Horse Radish Peroxidase (HRP).
特記事項 (標識)	Molar enzyme/ antibody protein ratio is 4:1
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab97155の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC		1/200 - 1/5000.
IHC-P		1/200 - 1/5000.
ELISA		1/10000 - 1/100000. ELISA-Primary
WB		1/5000 - 1/50000. WB-Colorimetric 1/5000 - 1/30000; WB-Chemiluminescent 1/10000 - 1/50000

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors