

Goat Anti-Human IgA alpha chain (HRP) preadsorbed ab98558

2 References

製品の概要		
製品名	Goat Anti-Human IgA alpha chain (HRP) preadsorbed	
由来種	Goat	
ターゲット生物種	Human	
アプリケーション	適用あり: ELISA, WB, IHC-P, ICC	
吸着処理血清	Mouse, Rat	
標識	HRP	
		more details

製品の特性	
製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Horse Radish Peroxidase.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション	
The Abpromise guarantee	<u>Abpromise保証は</u> 、次のテスト済みアプリケーションにおけるab98558の使用に適用されます アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		1/10000 - 1/100000. (as Primary).
WB		1/2000 - 1/25000. (Colorimetric: 1/2000 - 1/20000; Chemiluminescent: 1/5000 - 1/25000).

アプリケーション	Abreviews	特記事項
IHC-P		1/200 - 1/2000.
ICC		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors