

Goat Anti-Human IgE (SureLight® PE) ab99910

医薬用外毒物

製品の概要

製品名	Goat Anti-Human IgE (SureLight® PE)
由来種	Goat
ターゲット生物種	Human
アプリケーション	適用あり: Flow Cyt
標識	SureLight® PE

製品の特性

製品の状態	Lyophilized:rehydrate with 1.0 ml distilled deionized water
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at +4°C. Please see notes section. Store In the Dark.
バッファー	pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 1.64% Sodium phosphate, 1.7% Sucrose, 0.1% BSA, 0.87% Sodium chloride
精製度	Protein G purified
特記事項 (標識)	SureLight® R-phycoerythrin Excitation max: 565>498 nm Emission max.: 578 nm
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	Rehydrate with 1.0 ml distilled deionized water

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab99910の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use a concentration of 3 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors