

# Goat Anti-Human IgG H&L preadsorbed ab97166

## 2 References

### 製品の概要

---

|          |                                                                                                                                                                                                                                                                                                                         |
|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 製品名      | Goat Anti-Human IgG H&L preadsorbed                                                                                                                                                                                                                                                                                     |
| 由来種      | Goat                                                                                                                                                                                                                                                                                                                    |
| ターゲット生物種 | Human                                                                                                                                                                                                                                                                                                                   |
| 特異性      | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG and with light chains common to other Human immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. Reduced cross-reactivity to bovine, chicken, goat, horse, mouse, pig, rabbit and rat was detected. |
| アプリケーション | <b>適用あり:</b> ICC, IHC-P, ELISA, WB                                                                                                                                                                                                                                                                                      |
| 吸着処理血清   | Chicken, Cow, Horse, Mouse, Pig, Rabbit, Rat <a href="#">more details</a>                                                                                                                                                                                                                                               |

### 製品の特性

---

|          |                                                                                                                                                                                                                                         |
|----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 製品の状態    | Liquid                                                                                                                                                                                                                                  |
| 保存方法     | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.                                                                                                        |
| バッファー    | pH: 6.8<br>Preservative: 0.09% Sodium azide<br>Constituent: PBS                                                                                                                                                                         |
| 精製度      | Immunogen affinity purified                                                                                                                                                                                                             |
| 特記事項(精製) | Antiserum was cross adsorbed using bovine, chicken, horse, mouse, pig, rabbit and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. |
| ポリ/モノ    | ポリクローナル                                                                                                                                                                                                                                 |
| アイソタイプ   | IgG                                                                                                                                                                                                                                     |

### アプリケーション

---

The Abpromise guarantee

**Abpromise保証は、次のテスト済みアプリケーションにおけるab97166の使用に適用されます**

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション | Abreviews | 特記事項                                                                                     |
|----------|-----------|------------------------------------------------------------------------------------------|
| ICC      |           | Use at an assay dependent dilution.                                                      |
| IHC-P    |           | 1/200 - 1/2000.                                                                          |
| ELISA    |           | 1/100 - 1/10000. ELISA-Coating 1/100 - 1/500; ELISA-Primary 1/1000 - 1/10000             |
| WB       |           | 1/1000 - 1/30000. WB-Colorimetric 1/1000 - 1/10000; WB-Chemiluminescent 1/1000 - 1/30000 |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors