

Goat Anti-Human IgG+IgM+IgA H&L (HRP) ab102420

3 References

製品の概要

製品名	Goat Anti-Human IgG+IgM+IgA H&L (HRP)
由来種	Goat
ターゲット生物種	Human
特異性	By immunoelectrophoresis and ELISA this antibody reacts with human IgG, IgA and IgM. No antibody was detected against non-immunoglobulin serum proteins. This antibody may cross react with immunoglobulins from other species.
アプリケーション	<b>適用あり:</b> IHC-P, ELISA, WB, ICC
免疫原	Other Immunogen Type corresponding to IgG + IgM + IgA.
標識	HRP

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Constituents: 0.2% BSA, 0.05% CMIT/MIT based preservative
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to HRP.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee      Abpromise保証は、次のテスト済みアプリケーションにおけるab102420の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/200 - 1/5000.

アプリケーション	Abreviews	特記事項
ELISA		1/10000 - 1/100000. ELISA - Primary
WB		1/5000 - 1/50000. WB - Colorimetric: 1/5000 - 1/30000; WB - Chemiluminescent: 1/10000 - 1/50000
ICC		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors