

Goat Anti-Mouse IgG H&L (Alkaline Phosphatase) preadsorbed ab7069

★★★★☆ [1 Abreviews](#) [4 References](#) [画像数 1](#)

製品の概要

製品名	Goat Anti-Mouse IgG H&L (Alkaline Phosphatase) preadsorbed
由来種	Goat
ターゲット生物種	Mouse
アプリケーション	適用あり: ICC/IF, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, Dot blot, WB
吸着処理血清	Chicken, Cow, Goat, Guinea pig, Hamster, Horse, Human, Rabbit, Rat, Sheep more details
免疫原	Full length native mouse IgG (purified).
標識	Alkaline Phosphatase

製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.1% Sodium azide Constituents: 0.00136% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, 1% BSA
精製度	Affinity purified
特記事項 (精製)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using mouse IgG coupled to agarose beads.
特記事項 (標識)	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)
ポリ/モノ	ポリクローナル

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab7069の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/200 - 1/1000.
ELISA		1/7000.
IHC-P	★★★★★ (1)	Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.
Immunomicroscopy		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution.
WB		1/1000 - 1/2000.

画像



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Goat Anti-Mouse IgG H&L (Alkaline Phosphatase) preadsorbed (ab7069)
This image is courtesy of an anonymous Abreview.

ab13918 staining RANK in human skin tissue sections by Immunohistochemistry (IHC-P - paraformaldehyde-fixed, paraffin-embedded sections). Tissue was fixed with formaldehyde and blocked with 1% BSA/5% normal goat serum/0.1% saponin for 30 minutes at 22°C. Samples were incubated with primary antibody (1/50 in blocking buffer) for 9 hours at 4°C. An alkaline phosphatase-conjugated goat anti-mouse IgG H&L (ab7069) (1/300) was used as the secondary antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors