

Goat Anti-Mouse IgG2b heavy chain (PE/Cy7®) preadsorbed ab130790

6 References

製品の概要

製品名	Goat Anti-Mouse IgG2b heavy chain (PE/Cy7®) preadsorbed
由来種	Goat
ターゲット生物種	Mouse
特異性	ab130790 reacts with the heavy chain of Mouse IgG2b as demonstrated by ELISA, FLISA, and/or Flow Cytometry. Goat Anti-Mouse IgG2b heavy chain (PE/Cy7®) preadsorbed shows minimal cross reactivity with Human immunoglobulins. Cross absorption: Mouse IgM, IgG1, IgG2a, IgG3 and IgA, pooled Human sera and purified Human paraproteins.
アプリケーション	適用あり: ELISA, Flow Cyt, ICC/IF, WB
吸着処理血清	Human more details
免疫原	Full length protein corresponding to Mouse IgG2b. Mouse IgG2b paraproteins
標識	PE/Cy7®. Ex: 496nm, Em: 774nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	pH: 7.1 Preservative: 0.1% Sodium azide Constituents: PBS, Sucrose
精製度	Immunogen affinity purified
特記事項(精製)	ab130790 is purified by affinity chromatography on Mouse IgG2b covalently linked to agarose.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab130790の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration.
ICC/IF		Use 0.1µg for 10 ⁶ cells.
WB		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors