

Goat Anti-Rat IgG H&L (DyLight® 488) preadsorbed ab98420

5 References 画像数 2

製品の概要

製品名	Goat Anti-Rat IgG H&L (DyLight® 488) preadsorbed	
由来種	Goat	
ターゲット生物種	Rat	
アプリケーション	適用あり: Flow Cyt, ICC/IF, IHC-P	
吸着処理血清	Chicken, Cow, Goat, Horse, Human, Mouse, Rabbit, Sheep	more details
標識	DyLight® 488. Ex: 493nm, Em: 518nm	

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 488.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG


アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab98420の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		1/50 - 1/200.

アプリケーション	Abreviews	特記事項
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.


画像



Spectrum - Goat Anti-Rat IgG H&L (DyLight® 488)
preadsorbed (ab98420)

Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.



Immunocytochemistry/ Immunofluorescence - Goat
Anti-Rat IgG H&L (DyLight® 488) preadsorbed
(ab98420)

ICC/IF image of **ab6161** stained HeLa cells. The cells were 100% methanol fixed (5 min) and then incubated in 1%BSA / 10% normal goat serum / 0.3M glycine in 0.1% PBS-Tween for 1h to permeabilise the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody (**ab6161**, 1µg/ml) overnight at +4°C. The secondary antibody (green) was DyLight® 488 goat anti-rat- H&L, pre-adsorbed (ab98420) used at a 1/250 dilution for 1h. Alexa Fluor® 594 WGA was used to label plasma membranes (red) at a 1/200 dilution for 1h. DAPI was used to stain the cell nuclei (blue) at a concentration of 1.43µM.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors