

# Mouse monoclonal [B3102E8] Anti-Human IgE Fc (Biotin) ab99807

[7 References](#) [画像数 1](#)

### 製品の概要

製品名	Mouse monoclonal [B3102E8] Anti-Human IgE Fc (Biotin)
由来種	Mouse
ターゲット生物種	Human
特異性	ab99807 Reacts with the Fc portion of the heavy chain of Human IgE
アプリケーション	<b>適用あり:</b> ELISA, WB, Dot blot, IHC-Fr, ICC
標識	Biotin

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. Store In the Dark.
バッファー	pH: 7.4 Preservative: 0.1% Sodium azide Constituent: PBS
精製度	Affinity purified
ポリ/モノ	モノクローナル
クローン名	B3102E8
アイソタイプ	IgG1
軽鎖の種類	kappa
特記事項	This product has switched from ascites to tissue culture supernatant on 8 <sup>th</sup> June 2020. Lot numbers higher than GR3323141 are from tissue culture supernatant. Please note that the dilutions may need to be adjusted accordingly. If you have any questions, please do not hesitate to contact our scientific support team.

### アプリケーション


The Abpromise guarantee

Abpromise保証は、 次のテスト済みアプリケーションにおけるab99807の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		1/1000 - 1/2000.
WB		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.
ICC		Use at an assay dependent dilution.

## 画像



ELISA plate was coated with Human IgE. Immunoglobulin was detected with serially diluted Mouse monoclonal [B3102E8] Anti-Human IgE Fc (Biotin) (ab99807) followed by Streptavidin-HRP and ABTS.

This product has switched from ascites to tissue culture supernatant on 8<sup>th</sup> June 2020. Lot numbers higher than GR3258714 are from tissue culture supernatant.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors