

## Product datasheet

# Mouse monoclonal [B3506B4] Anti-Human IgA1 Fc (PE) ab99797

### 製品の概要

---

製品名	Mouse monoclonal [B3506B4] Anti-Human IgA1 Fc (PE)
由来種	Mouse
ターゲット生物種	Human
特異性	ab99797 reacts with the Fc portion of the heavy chain of IgA1.
アプリケーション	<b>適用あり:</b> ELISA, Dot blot, WB, IHC-Fr, ICC/IF
標識	PE. Ex: 488nm, Em: 575nm

### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	pH: 7.3 Preservative: 0.1% Sodium azide Constituents: Sucrose, PBS
精製度	Protein G purified
ポリモノ	モノクローナル
クローン名	B3506B4
アイソタイプ	IgG1
軽鎖の種類	kappa
特記事項	Source: Ascites fluid

### アプリケーション

---

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab99797の使用に適用されず**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution. for Dot and slot immunoblotting
WB		Use at an assay dependent dilution. Predicted molecular weight: 38 kDa.
IHC-Fr		Use at an assay dependent dilution.
ICC/IF		Use 0.2µg for 10 <sup>6</sup> cells.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors