

### Rabbit Anti-Goat IgG Fc (HRP) preadsorbed ab102345

#### 製品の概要

製品名	Rabbit Anti-Goat IgG Fc (HRP) preadsorbed	
由来種	Rabbit	
ターゲット生物種	Goat	
アプリケーション	適用あり: ICC, IHC-P, ELISA, WB	
吸着処理血清	Human, Mouse, Rat	<a href="#">more details</a>
標識	HRP	

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Constituent: 0.05% CMIT/MIT based preservative
精製度	Immunogen affinity purified
特記事項 (精製)	Antiserum was solid phase adsorbed to ensure class specificity. Antiserum was cross adsorbed using human, mouse and rat immunosorbents to remove cross reactive antibodies. The antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Horse Radish Peroxidase.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	By immunoelectrophoresis and ELISA this antibody reacts specifically with goat IgG. Cross reactivity with IgA and IgM is negligible. No antibody was detected against non-immunoglobulin serum proteins. Less than 1% cross reactivity to human, mouse and rat IgG was detected. This antibody may cross react with IgG from other species.

#### アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab102345の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC		1/200 - 1/5000.
IHC-P		1/200 - 1/5000.
ELISA		1/10000 - 1/100000. (Primary).
WB		1/5000 - 1/50000. (Colorimetric: 1/5000 - 1/30000; Chemiluminescent: 1/10000 - 1/50000).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors