

### Rabbit Anti-Human IgA alpha chain (DyLight® 550) ab97000

**3 References** [画像数 1](#)

#### 製品の概要

製品名	Rabbit Anti-Human IgA alpha chain (DyLight® 550)
由来種	Rabbit
ターゲット生物種	Human
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgA. Cross reactivity with other immunoglobulins and light chains is less than 0.1%.
アプリケーション	<b>適用あり:</b> IHC-P, ICC/IF, Flow Cyt
標識	DyLight® 550. Ex: 562nm, Em: 576nm

#### 製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS
精製度	Immunogen affinity purified
特記事項 (精製)	Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 550.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	DyLight® 550 replaces DyLight® 549.

#### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、**次のテスト済みアプリケーションにおけるab97000の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/500. DyLight® 550 replaces DyLight® 549. The DyLight® 550 excitation is 562nm and the emission is 576nm. DyLight® 550 provides the same orange-to-red fluorescence and photostability as the DyLight® 549 dye. Please refer to the vial in order to check whether the product you received is DyLight® 549 or DyLight® 550. This information however will have no impact on your
ICC/IF		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

## 画像



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors