

### Rabbit Anti-Human IgM mu chain (Biotin) ab97208

#### 製品の概要

|          |  |
|----------|--|
| 製品名      | Rabbit Anti-Human IgM mu chain (Biotin)  |
| 由来種      | Rabbit   |
| ターゲット生物種 | Human  |
| 特異性      | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgM. Cross reactivity with other immunoglobulins and light chains is less than 0.1%. |
| アプリケーション | <b>適用あり:</b> IHC-P, ICC/IF, ELISA, WB  |
| 標識       | Biotin   |

#### 製品の特性

|           |   |
|-----------|---|
| 製品の状態     | Liquid  |
| 保存方法      | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.   |
| バッファー     | pH: 6.8<br>Preservative: 0.09% Sodium azide<br>Constituents: 0.2% BSA, PBS  |
| 精製度       | Immunogen affinity purified   |
| 特記事項 (精製) | Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Biotin. Biotinylated antibody was demonstrated with avidin/peroxidase when coated on microtiter wells. |
| ポリ/モノ     | ポリクローナル   |
| アイソタイプ    | IgG   |

#### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab97208の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項            |
|----------|-----------|-----------------|
| IHC-P    |           | 1/250 - 1/2500. |

| アプリケーション | Abreviews | 特記事項                                    |
|----------|-----------|---|
| ICC/IF   |           | 1/100 - 1/500.                          |
| ELISA    |           | 1/10000 - 1/200000. ELISA-Primary       |
| WB       |           | 1/10000 - 1/200000. WB-Chemiluminescent |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors