

Rabbit Anti-Mouse IgG H&L (Biotin) preadsorbed ab98784

★★★★★ [1 Abreviews](#) [画像数 1](#)

製品の概要

製品名	Rabbit Anti-Mouse IgG H&L (Biotin) preadsorbed	
由来種	Rabbit	
ターゲット生物種	Mouse	
アプリケーション	適用あり: IHC-P, ICC/IF, ELISA, WB	
吸着処理血清	Human, Rat	more details
標識	Biotin	

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Biotin.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG


アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab98784の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/500 - 1/5000.

アプリケーション	Abreviews	特記事項
ICC/IF		1/200 - 1/1000.
ELISA		1/20000 - 1/400000. as primary
WB		1/20000 - 1/400000. chemiluminescent detection

画像



Immunohistochemical analysis of dog skin tissue, detecting staining of NQO1 with ab98784.

Tissue was fixed with formaldehyde and blocked with none for ; antigen retrieval was by heat mediation in a citrate buffer (pH 6). Samples were incubated with primary antibody (**ab28947 Mouse monoclonal [A180] to NQO1** at 1/175 in BSA in TBS) for 45 minutes. ab98784 (1/500) was used as the secondary antibody.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Rabbit Anti-Mouse IgG H&L (Biotin) preadsorbed (ab98784)

This image is courtesy of an anonymous Abreview

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors