

Rabbit Anti-Rat IgG H&L (beta-galactosidase) ab136716

製品の概要

製品名	Rabbit Anti-Rat IgG H&L (beta-galactosidase)
由来種	Rabbit
ターゲット生物種	Rat
アプリケーション	<b>適用あり:</b> IHC-Fr, IHC-P, WB, ELISA, ICC
免疫原	Highly purified whole immunoglobulin from Rat.
標識	beta-galactosidase

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Constituents: 49% PBS, 50% Glycerol
精製度	Immunogen affinity purified
特記事項 (精製)	ab136716 is purified by a very stringent purification protocol where non-specific and specific purifications alternate (affinity chromatography).
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	Preparation of Substrate Solutions: 0.1M phosphate buffer pH7, 1mM MgSO <sub>4</sub> , 0.2mM MnSO <sub>4</sub> , 2mM Mg titriplex. Add extemporaneously 7µl/ml beta-mercaptoethanol and 0.9mg/ml ortho-nitrophenyl-beta-galactopyranoside. Incubate for 30 minutes at 37°C. Stop the reaction with V/2 volume of 1M Na <sub>2</sub> CO <sub>3</sub> .

アプリケーション

The Abpromise guarantee      **Abpromise保証は、次のテスト済みアプリケーションにおけるab136716の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-Fr		Use at an assay dependent concentration.

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.
ELISA		1/100 - 1/1000.
ICC		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors