

Rabbit Anti-Sheep IgG H&L (Glucose Oxidase) ab136721

製品の概要

製品名	Rabbit Anti-Sheep IgG H&L (Glucose Oxidase)
由来種	Rabbit
ターゲット生物種	Sheep
アプリケーション	適用あり: IHC-P, IHC-Fr, ELISA, ICC, WB
免疫原	Highly purified whole immunoglobulin from Sheep.
標識	Glucose Oxidase

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Constituents: 49% PBS, 50% Glycerol
精製度	Immunogen affinity purified
特記事項 (精製)	ab136721 is purified by a very stringent purification protocol where non-specific and specific purifications alternate (affinity chromatography).
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	Do not store in plastic tubes (i.e. eppendorf) to avoid losing activity. Preparation of Substrate Solutions: ELISA: 10mg/ml D(+)-glucose, 0.1mg/ml grade II (Rz=0.6) peroxidase, 0.35mg/ml OPD in 0.1M pH 5 citrate-phosphate buffer. Incubation 30 minutes at 37°C. Stop the reaction with 1/2 volume of 2N H ₂ SO ₄ . Reading at 492 nm. Immunocytochemistry: 15mg/ml D(+)-glucose, 0.5mg/ml blue thiazolyl, 0.35mg/ml phenazine methosulfonate to be dissolved in that order in a 0.1M pH 6.8 phosphate buffer. Incubation 20 to 30 minutes in darkness.

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab136721の使用に適用されず
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/10 - 1/100.
IHC-Fr		1/10 - 1/100.
ELISA		1/100 - 1/1000.
ICC		1/10 - 1/100.
WB		1/400.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors