

Rabbit Anti-Sheep IgG H&L (HRP) preadsorbed ab7111

5 References

製品の概要

製品名	Rabbit Anti-Sheep IgG H&L (HRP) preadsorbed
由来種	Rabbit
ターゲット生物種	Sheep
アプリケーション	適用あり: ICC/IF, IHC-P, WB, ELISA, Immunomicroscopy, Dot blot, IHC-Fr
吸着処理血清	Human more details
免疫原	Full length native Sheep IgG (purified).
標識	HRP

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.50 Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride, 1% BSA
精製度	Affinity purified
特記事項(精製)	Horseradish Peroxidase (HRP)
特記事項(標識)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
ポリ/モノ	ポリクローナル

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab7111の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/300 - 1/1500.
IHC-P		Use at an assay dependent dilution.
WB		1/500 - 1/3000.
ELISA		1/5000 - 1/22000.
Immunomicroscopy		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors