

## Product datasheet

### DU 145 nuclear extract lysate ab14645

#### 製品の概要

---

製品名	DU 145 nuclear extract lysate
特記事項	We recommend aliquoting the extracts into single-use fractions and then storing them.

#### 製品の特性

---

Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.90 Constituents: 0.75% Potassium chloride, 0.009% PMSF, 0.01% Magnesium chloride, 0.008% DTT, 0.48% HEPES, 20% Glycerol
背景	DU 145 cells are used as an in vitro model for prostate cancer. They grow in continuous culture as a monolayer, are used for inducing adenocarcinoma in nude mice and are weakly positive for acid phosphatase.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors