

Rat brain tissue lysate - mitochondrial extract ab110348

製品の概要

製品名	Rat brain tissue lysate - mitochondrial extract
特記事項	<p>Product was previously marketed under the MitoSciences sub-brand.</p> <p>2 mg of rat brain mitochondria in 400ul of iso-osmotic buffer. This sample can act as a positive control in immunocapture applications and Blue Native PAGE.</p> <p>1 mg solubilized mitochondria for Immunocapture experiments. For the Immunocapture of OXPHOS complexes it is recommended that 180 µl of the mitochondrial membranes supplied be solubilized by the addition of 20 µl of 10% lauryl maltoside. Insoluble material is removed by centrifugation at 72,000g before immunocapture of complexes from the supernatant.</p>
アプリケーション	適用あり: IP, SDS-PAGE, WB

製品の特性

Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
バッファー	<p>pH: 7.5</p> <p>Constituents: 0.02% PMSF, 0.16% Tris, 0.89% EDTA, 8.56% Sucrose</p>
ライゼート 備考	<p>2 mg of rat brain mitochondria in 400ul of iso-osmotic buffer. This sample can act as a positive control in immunocapture applications and Blue Native PAGE.</p> <p>1 mg solubilized mitochondria for Immunocapture experiments. For the Immunocapture of OXPHOS complexes it is recommended that 180 µl of the mitochondrial membranes supplied be solubilized by the addition of 20 µl of 10% lauryl maltoside. Insoluble material is removed by centrifugation at 72,000g before immunocapture of complexes from the supernatant.</p>

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab110348の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IP		Use at an assay dependent concentration.
SDS-PAGE		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors