

## Product datasheet

# Rabbit Anti-Chicken IgY H&L (Biotin) ab6752

9 References [画像数 2](#)

### 製品の概要

製品名	Rabbit Anti-Chicken IgY H&L (Biotin)
由来種	Rabbit
ターゲット生物種	Chicken
アプリケーション	<b>適用あり:</b> IHC-Fr, Immunomicroscopy, WB, Dot blot, ELISA, IHC-P, ICC/IF
免疫原	chicken IgY whole molecule
標識	Biotin

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium azide Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Potassium phosphate
精製度	Affinity purified
特記事項 (精製)	Biotinamidocaproate N-Hydroxysuccinimide Ester (BAC). Biotin/Protein Ratio: 10-20 BAC molecules per Rabbit IgG molecule.
特記事項 (標識)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using chicken IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
ポリ/モノ	ポリクローナル


### アプリケーション

**The Abpromise guarantee** [Abpromise保証は](#)、次のテスト済みアプリケーションにおけるab6752の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-Fr		Use at an assay dependent concentration.


アプリケーション	Abreviews	特記事項
Immunomicroscopy		Use at an assay dependent concentration.
WB		1/2000 - 1/20000.
Dot blot		Use at an assay dependent concentration.
ELISA		1/500000.
IHC-P		Use at an assay dependent concentration.
ICC/IF		1/1000 - 1/5000.

## 画像



Immunocytochemistry/ Immunofluorescence - Rabbit  
Anti-Chicken IgY H&L (Biotin) (ab6752)

Ab6752 was used with the primary antibody [ab13989](#) in ICC/IF.  
See Abreview on [ab13989](#).



Immunohistochemistry (Frozen sections) - Rabbit  
Anti-Chicken IgY H&L (Biotin) (ab6752)

Ab6752 was used at dilution 1/1000 with the primary antibody  
[ab13970](#) in IHC-Fr. See the review on [ab13970](#).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors