

# Recombinant Human Aprataxin protein ab93630

画像数 1

### 製品の詳細

製品名	Recombinant Human Aprataxin protein
精製度	> 95 % SDS-PAGE. ab93630 is purified using conventional chromatography techniques.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWAGSMQD PKMQVYKDEQ VVVIKDKYPK ARYHWLVLPW TSISSLKAVA REHLELLKHM HTVGEKVIVD FAGSSKLRFR LGYHAIPSMS HVHLHVISQD FDSPCLKNKK HWNSFNTEYF LESQAVIEMV QEAGRVTVRD GPELLKLPL RCHECQQLLP SIPQLKEHLR KHWTQ
領域	1 to 168

### 特性

Our **Abpromise guarantee** covers the use of **ab93630** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid


### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  pH: 7.50 Constituents: 0.00174% PMSF, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
------------	---

## 関連情報

機能	DNA-binding protein involved in single-strand DNA break repair, double-strand DNA break repair and base excision repair. Resolves abortive DNA ligation intermediates formed either at base excision sites, or when DNA ligases attempt to repair non-ligatable breaks induced by reactive oxygen species. Catalyzes the release of adenylate groups covalently linked to 5'-phosphate termini, resulting in the production of 5'-phosphate termini that can be efficiently rejoined. Also able to hydrolyze adenosine 5'-monophosphoramidate (AMP-NH(2)) and diadenosine tetraphosphate (AppppA), but with lower catalytic activity.
組織特異性	Widely expressed. In brain, it is expressed in the posterior cortex, cerebellum, hippocampus and olfactory bulb. Isoform 1 is highly expressed in the cerebral cortex and cerebellum, compared to isoform 2.
関連疾患	Defects in APTX are the cause of ataxia-oculomotor apraxia syndrome (AOA) [MIM:208920]. AOA is an autosomal recessive syndrome characterized by early-onset cerebellar ataxia, oculomotor apraxia, early areflexia and late peripheral neuropathy. Defects in APTX are a cause of coenzyme Q10 deficiency (COQ10D) [MIM:607426]. Coenzyme Q10 deficiency is an autosomal recessive disorder with variable manifestations. It can be associated with three main clinical phenotypes: a predominantly myopathic form with central nervous system involvement, an infantile encephalomyopathy with renal dysfunction and an ataxic form with cerebellar atrophy.
配列類似性	Contains 1 C2H2-type zinc finger. Contains 1 FHA-like domain. Contains 1 HIT domain.
ドメイン	The histidine triad, also called HIT motif, forms part of the binding loop for the alpha-phosphate of purine mononucleotide. The FHA-like domain mediates interaction with NCL; XRCC1 and XRCC4. The HIT domain is required for enzymatic activity. The C2H2-type zinc finger mediates DNA-binding.
細胞内局在	Nucleus > nucleoplasm. Nucleus > nucleolus. Upon genotoxic stress, colocalizes with XRCC1 at sites of DNA damage. Colocalizes with MDC1 at sites of DNA double-strand breaks. Interaction with NCL is required for nucleolar localization.

## 画像



15% SDS-PAGE showing ab93630 at approximately 23.9kDa (3µg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors